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© all rights reserved, Qingdao Chuangmeng Instrument Co., Ltd

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Please read the Instruction Manual carefully, for correctly grasping the installation and using method of this product. Please keep properly this Instruction Manual after reading, for the usage during troubleshooting and maintenance in the future.

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Address: No. 7 Wenzhou Road, City Northern District, Qingdao City, China

Production base: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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1. Summary

Torsional spring dynamometer is a special testing instrument designed to correct the stiffness and indication error of torsion spring in rotary viscometer. It is also an effective instrument for regularly measuring torsion spring to ensure the accuracy of viscometer.

2. Model and specification

Model	Name	Configuration
1105	Torsional spring dynamometer	

3. The main technical parameters

Name	Technical parameter
Measuring range	Within the range of standard value of the rigidity of the spring, can determine in succession (Weight: 1g—70g)
External size	200×142×300mm

4. The Structure and Working Principles of the Instrument

The structure is mainly composed of the following parts

1) Base: composed of base and double vertical shaft. The height of the support is regulable and located. They form the main body and can make dynamometer align viscometer corresponding position.


2) Runner part: it is used to connect the inner cylinder and the weight as the guide part to determine the rotation angle of the torsion spring on the viscometer under the measured load.

3) working principle

1. When measuring with torsion spring dynamometer, the measuring element inner cylinder of six speed viscometer is used as the force transmission medium, and the correction method is direct.

2. The height of the viscometer is adjusted automatically by the support, and the height of the inner cylinder of the viscometer is aligned, and the positioning is reliable.

3. Working diagram of dynamometer:

(Figure 2) Working diagram of dynamometer

Number	Part number	Name and specification	Quantity	Remarks
1	1105	Viscometer torsion spring	1	
2		Viscometer inner cylinder shaft	1	
3		Dynamometer pulley holder	1	
4		Traction line	1	With hook
5		Dynamometer pulley	1	
6		Weight	1	
7		Viscometer inner cylinder	1	

(Table 2) Working schematic structure list

5. Operation of apparatus

- 1.Remove the instrument to be tested (viscometer) and the outer rotating cylinder. Level with dynamometer.
- 2.Make the reel drum with the inner of the viscometer. Fix one end of the draw wire with the adhesive tape to the middle of the inner, circle the inner towards the left then draw straight and circle the pulley on the support .Hung up the weights on the other end.
- 3.Loosen the adjusting knob (10) (see fig.1) and adjust the height of the axis (4) in order to make the height of the rotary wheel of the dynamometer alignment to the inner of the viscometer. Then begin to do the efficacy of the torsional spring.
- 4.Adjust the weights according to “rigidity linear test form of torsional spring “.Read the measured value from the viewing window and check the torsional spring.

Computational method of measuring standard value,

$$KS = \frac{G.r}{\phi}$$

Gr---weights (g)


r---radius of inner (r=1.7245)

φ---dial reading (degree)

K_S---Standard spring rigidity

Example: A R₁-B₁ measuring group. Under standard spring rigidity value (F1 torsion spring:K_S=0.394g. centimeter per degree), add the weights 50g then the corner of the needle is:

$$\phi = \frac{50 \times 1.7245}{0.394} = 218.8 \text{ (degree)}$$

(Figure 3) Operation schematic diagram

Measures with different weights should conform to the torsion value shown in the ‘rigidity linear test form of spring’

Weights (g)	Scope (case)	Actual measurement (case)	Weights(g)	Scope (case)	Actual measurement (case)
5	21.55~22.21		40	172.45~177.70	
10	43.11~44.42		45	194.00~199.90	
15	64.66~66.63		50	215.57~222.13	
20	86.23~88.85		55	237.12~244.30	
25	107.78~111.06		60	258.68~266.56	
30	129.34~133.28		65	280.23~288.77	
35	150.89~155.49		70	301.80~310.99	

(Table 3) F1Rigidity linear test form of torsional spring

Weights (g)	Scope (case)	Actual measurement (case)	Weights(g)	Scope (case)	Actual measurement (case)
1	21.55~22.21		8	172.45~177.70	
2	43.11~44.42		9	194.00~199.90	
3	64.66~66.63		10	215.57~222.13	
4	86.23~88.85		11	237.12~244.30	
5	107.78~111.06		12	258.68~266.56	
6	129.34~133.28		13	280.23~288.77	
7	150.89~155.49		14	301.80~310.99	

(Table 4) Rigidity linear test form of torsional spring

6. The maintenance of the instrument

- 1). The operator should be familiar with all the operation process and possibly occurring circumstances, and operate according to the testing procedures
- 2). Having finished operating, each part should be clean out. Take down the weights, put and keep properly alone
- 3). Take and place carefully when moving, repairing and cleaning the instrument in order to avoid causing parts deformed and influencing precision and use
- 4). As fig.1 shows, between axis (4) and sleeve (2), pulley (8) and fixed axis, lubricating oil must be applied frequently.

7. The Transportation and Storage of Instrument

The transportation and storage of instrument should correspond to the JB/T9329-1999 standard. The product should be stored in the room with ventilation. The indoor air does not contain the impurity which can arouse device corrosion.

8. Fault judge and remove

Sequence number	Fault	Cause	Removing method
1	As it shown in Fig.1 the pulley (8) can not rotate flexibly.	The friction drag is produced between pulley (8) and fixed axis	Lubricating oil should be applied as required between pulley(8) and fixed axis
2	As it shown in Fig.1,the axis(4)and sleeve(2) can not rotate flexibly	The friction drag is produced between pulley (4) and sleeve(4)	Lubricating oil should be applied as required between axis(8) and sleeve(2)

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Production address:No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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Model of the main motor:1105

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Manufacturing No:

No	Part number	Name	Unit	Quantity	Remarks
1	P03A2	砝码（100g） Gweight（100g）	set	1	
2	P03A3	牵引线 Traction	piece	1	The line is 1 meter long
3		使用说明书 Manual		1	

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Qingdao Chuangmeng Instrument Co., Ltd. Packing list
Manufacturing enterprise: Qingdao Chuangmeng Instrument Co.,Ltd.