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Qingdao Chuangmeng Instrument Co., Ltd

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Frequency Conversion HighSpeed Mixers

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© all rights reserved, Qingdao Chuangmeng Instrument Co., Ltd

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Carefully read this User Manual to learn how to install and use the product correctly. After reading, properly keep the User Manual as a reference for future maintenance and repair.

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Address: No. 7 Wenzhou Road, City Northern District, Qingdao City, China

Production base: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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I .Summary

The frequency conversion high speed mixer adopts the advanced frequency conversion technology and microcomputer technology at present. It has functions of timing, constant speed, single shaft and multi axis simultaneous operation. Mixing blades were made by API standard. It has the characteristics of simple operation, large starting torque, stable speed, reliability and low noise. It is a special experimental equipment for petroleum, exploration and chemical industry. Mainly used for high-speed stirring various oil drilling fluids, chemical reagents, etc., is an ideal high-speed stirring equipment.

II.Type and specification

Model	Parts
1613	1 electric motor ,1 liquid cup.
1614	2 electric motor ,2 liquid cup 。
1615	4 electric motor ,4 liquid cup.

III.Parameters of the instrument


Parts Model Name	1613	1614	1615
Working power supply	220V±5% AC 50/60Hz		
Rated power	200 W	200W×2	200W×4
Unloaded speed	3000、4000、6000、8000、10000、11000、12000r/min		
Stirring volume	350ml	350ml×2	350ml×4
Timing range	1~99min		
Atemperature	~40℃		
Ambient humidity	10~85RH%		

IV. Structure and working principle


(1)Mixer.

(2) Control box


1.Schematic diagram of Model 1613 high speed mixer control box


2.Schematic diagram of Model 1614 high speed mixer controlbox

3. Schematic diagram of Model 1615 high speed mixer control box

(3) Structure chart and schedule

No	Name	No	Name
1	Shaft as-mixing	3	Housing
2	Liquid cup	4	Control box

V. Instrument operation

1. Connect the instrument to the corresponding voltage and power supply.

Note: the four axis mixer has two output lines, which should be connected to their respective mainframe respectively.

2. Pour the sample into the liquid cup. In order to prevent overflow when sample is stirred, the injection amount of sample is generally not more than 2/3.

3. Place the top of the liquid cup on the opening of the cup assembly, and place the cup on the lower support seat. The liquid cup is held tightly by the pressing cup assembly and the lower supporting seat.

4. When selecting the operation, select the gear according to the tachometer and set the mixing time. Finally, press the gear key to start the operation of the instrument, and the gear number flashes during the operation.

5. Select stop, (a) press the gear knob. (b) set the time until the instrument stops automatically. Note: gear position and timing can only be adjusted when the instrument is stopped.

6. After stirring, hold the cup tightly and lift it forcefully. Remove the cup from the lower support.

Comparison table between gear and speed:

Gear	1	2	3	4	5	6	7
Speed	3000	4000	6000	8000	10000	110000	120000

VI.Maintenance and maintenance.

1. Clean the components and dry them up. The instruments are placed in a dry environment.
2. When moving, repairing or maintaining instruments. It is necessary to take lightly and put it lightly so as not to cause deformation of parts and affect accuracy and use.
- 3.The stirring shaft must not collide, so as to avoid bending deformation and damage the instrument.
- 4.Motor and speed regulator should be kept clean, dry and not damp. Regular inspection should be carried out. It is strictly forbidden to store corrosive chemicals with the machine.
- 5.If the weight of the stirring blade is less than 5.1g, it needs to be replaced. Its original weight is about 5.5g.

VII.Judgement and elimination of faults

Fault	Reason	Maintenance method
Motor noise increase	Motor damage	Replace motor
Motor does not rotate	1. Motor damage 2. Circuit board failure.	1. Replace the motor 2. Replacing circuit boards
Only partial rotation speed	"Speed conversion button" damage	Change the speed conversion button.
Stirring blades do not turn.	Loosening of fastening screws	Tighten fastening screws.

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Production address:No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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Model of the main motor:Frequency Conversion High Speed Mixer

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Qingdao Chuangmeng Instrument Co., Ltd. Packing list
Manufacturing enterprise: Qingdao Chuangmeng Instrument Co.,Ltd.